

DATA SHEET E 260

DIAPHRAGM VACUUM PUMP LABOPORT[®] N 810.3 FT.18



Technical features

- 100% oil-free transfer thereby pure transferring, evacuation and compression of gases
- High level of vapor and condensate compatibility
- TFM[™] PTFE pump head and PTFE-coated diaphragm are ideal for extremely aggressive/corrosive gases and vapors
- (Ex) ATEX compliant according to
 ATEX II 2G IIB+H2 T3X internal atmosphere only
- Environmentally friendly
- High level of gas tightness
- Also available with an integrated gas ballast valve: N 810.3 FT.18 G

Technical data

Flow rate [m ³ /h] at atm. pressure	0.6	
Flow rate [l/min] at atm. pressure	10	
Ultimate vacuum [mbar abs.]	8.0 ≤ 30.0 with gas ballast valve in use	
Operating pressure [bar]	1	
Voltage [V]	230	
Frequencies [Hz]	50	
Motor protection	IP 44	
Power P ₁ [W]	90	
Operating current [A]	0.6	
With thermal switch and power fuse		
Motors with other voltages and frequencies on request		

MATERIALS

Pump head	TFM™ PTFE
Diaphragm	PTFE-coated
Valves	FFPM


OTHER DATA

Weight [kg]	6.9
Hose connections	ID 10
Dimensions W x H x D [mm]	140 x 187 x 281
Permissible media and ambient temperature [°C]	+5 to +40

ACCESSORIES	ID
Spares kit	057357



Dimensions



Performance characteristics

